



ES-SA002DE13

Dental Stool - Doctor Chair



Description:

- No wooden material inside the seat. The seating structure is made by plastic injection molding as a whole unit.
- No screw hole on the backrest.
- Conformity to the human ergonomics.
- Guarantee: 1 Year

Specification:

- Size: 45 x 52 x 84-98 cm (W x D x H)
- Upholstery / Leather: PU foaming /PVC
- Wheels: Translucence PU
- Base: Aluminum

Function:

- Height Adjustable
- Backrest Adjustable
- Seat Tiltable

Optional:

- Hand Stitch Micro fiber Leather
- Higher distance Gas Lift
- Footring



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com





ES-SA002DE13

Tabouret dentaire - Chaige pour Docteur



Description:

- Pas de matériel en bois dans le siège. La structure de l'assise est faite par moulage à injection de matière plastique comme unité complète.
- Pas de vis dans le dossier.
- Conforme à l'ergonomie de l'homme.
- Garantie: 1 an

Spécification:

- Taille: 45 x 52 x 84-98 cm (L x P x H)
- Rembourrage / Cuir: mousse PU / PVC
- Roues: PU Translucen e
- Base: Aluminium

Fonction:

- H
-)
- S

9b`Option:

- Cuir cousu à main
- Distance supérieure du vérin à gas
- Repose-pieds



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

