



# ESA H+\$

Precision Rotary  
Microtome



LABORATORY



# ESAH+\$

## Precision Rotary Microtome

It is modern and reliable sectioning device designed for precise sectioning of tissues up to 1 micron thickness.

The feed indicator, located in front of microtome is operated by the 'CAM DRIVE SYSTEM' for accurate feed setting from 1-50 microns in steps of one micron.

The interior mechanism rests on a heavy cast iron base and is covered with a full swing protective cover for easier cleaning and lubrication.

The salient features include independent feed mechanism with automatic safety device, universal knife holder with lateral movements permitting the use of the entire knife edge, universal vice type object holder for accurate centring of the specimen.

The instrument ensures life-long smooth operation.



**Cod. ESMT70**

### Technical Features

- Reversing Handle
- Pull-up feed
- Scale Knob
- Universal vice
- Universal Razor Holder
- Digital Scale
- Slides on Bali bearing

### Standard Accessories

- 1 x Razor 120 mm, with back & handle in box
- 1 x Oil Can
- 1 x Honing Stone in box
- 1 x Dust Cover
- 1 x Object holders (set of three)
- 1 x Disposable Blade Holder with 5 disposable Blades
- 1 x Illustrated operating manual



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

