

MICROPIPETTE OP LINE



The pipettes use principal of suction and drainage, control moving distance of the piston in its sleeve, negative pressure occurred, the quantitative liquid is drew.

The moving distance of piston is controlled by adjustable wheel. It is used conveniently and reliably. Especially the pipettes under volume 200 μ l as the tips used are the same, so it avoids the trouble to purchase many types of tip.

The modern line of the Pipetter series OP LINE 7 models for volumes from 0.2 to 5000 μ l.

The Function buttons for the pipetting, tip ejection and volume setting are separated to prevent errors. The clear and Large display shows the setted volume.

MODEL	VOLUME	% TOLERANCE	REPEATABILITY
OP2	0.2	12.0	<6.00
	0.5	5	<2.50
	2	1.5	<0.70
OP10	1	2.5	<1.25
	5	1.5	<0.60
	10	1.0	<0.40
OP20	2	5.0	<1.50
	10	1.0	<0.50
	20	1.0	<0.30
OP100	10	2.0	<0.60
	50	0.8	<0.24
	100	0.8	<0.15
OP200	20	1.2	<0.50
	100	0.8	<0.25
	200	0.8	<0.15
OP1000	100	1.8	<0.40
	500	0.8	<0.25
	1000	0.8	<0.15
OP5000	1000	1.2	<0.30
	2000	0.6	<0.25
	5000	0.6	<0.16



ESSE 3 Via Garibaldi 30
14022Castelnuovo D.B. (AT)
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.4ci.com





FIXED AND VARIABLE MICROPIPETTES

