

ESF5-100 SERIE SURGICAL ASPIRATORS

CE
0123



ESF5- 0 SERIE SURGICAL ASPIRATORS

ASPIRATEURS DE CHIRURGIE

Mobile suction machine for professional use, with powerful double piston motor.
ABS anti-shock case.
Four castors, two of which with brakes.
High vacuum 750 mm/Hg.
Vacuum gauge for the control of vacuum.
Vacuum regulation.
Two jars, 2 Lt. capacity made of Makrolon autoclavable.
Each jar with double security valve.
Indipendente double receiver jar.
Filter disposable.

TECHNICAL SPECIFICATIONS

- Suction power : 100 lt./min.
- Maximum vacuum : 750 mm/hg.
- Power consumption: 270 VA
- Two safety fuse
- Voltage: 230V - 50/60 Hz
- Size: 480x460x420(H) mm
- Weight: 18,0 kg

Aspirateur mobile à usage professionnel avec moteur puissant a double piston.
Boîtier en ABS antichoc.
Équipé avec quatre roulettes, dont 2 avec freins.
Vide maximum 750 mm/Hg.
Manomètre de contrôle de vide.
Régulateur de dépression.
Deux bocal 2 Lt. en Makrolon, autoclavables.
Chaque bocal est équipé d'une double valve de sécurité.
Fonctionnement autonome de chaque bocal.
Filtre usage unique,

CARACTÉRISTIQUES TECHNIQUES

- Puissance d'aspiration : 100 lt./min.
- Aspiration maximum: 750 mm/hg.
- Puissance absorbée: 270 VA
- Deux fusibles de sûreté.
- Alimentation: 230V - 50/60 Hz
- Dimensions: 480x460x420(H) mm
- Poids: 18,0 kg

SURGICAL TRAY



Available with :

- Double jar 2x2 Lt
- Double jar 2x4 Lt

Optional accessories

- Antistatic wheels
- Pneumatic footswitch
- Suction curette

Disponibile avec :

- Double bouteille 2x4 Lt
- Double bouteille 2x4 Lt

Accessoires optionnels

- roues antistatiques
- Pédale pneumatique
- Sonde aspiration



ESSE 3 Via Garibaldi 30
14022Castelnuovo D.B. (AT)
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.esse3-medical.com

