



# DATA SHEET

## Model: CP-400

# Ophthalmic chart projector

50 charts and various masks  
 Landolt, Letter, Child, Snellen, 50 kinds of charts and horizontal line, vertical line, one chart, various masks can be selected by the remote control box  
 Fast, silent chart disk rotation.  
 It moves to the following chart in a maximum of 0.3 seconds.

User program mode:  
 2kinds of programs can be selected 48 kinds(maximum) of chart and mask are memorisable.

1. Projection Distance: 2.0 ~ 7.0m
2. Projection magnification: 30  
(when the projection distance is 5m)
3. Projection size: 335mm X222mm  
(when the projection distance is 5m)
4. Tilt angle: from the horizon up or down for each 10 degrees
5. source: 12V 35W (halogen)
6. Dimensions: 248x336x234 mm (WxLxH)
7. Weight: 4kgs
8. Power Supply 100-240V
9. 120VA 5 minutes after the power automatically turn off

Standard:  
 remote control board, SUM4 alkaline batteries 2, screen, test visual acuity drawings, dust cover, spare bulb, spare fuse

Optional:  
 Replacement Parts: table sets, wall mounting brackets, visual acuity chart projector bracket, polarization color filter head, red and green lens, polarizer






ESSE3 srl, Via Garibaldi 30  
 14022 Castelnuovo D.B. (AT)  
 Tel +39 011 99 27 706  
 Fax +39 011 99 27 506  
 e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
 web: [www.esse3-medical.com](http://www.esse3-medical.com)

